



BME

BRIBIE MEDICAL EQUIPMENT

*Local Family
Operated Business*

PRODUCT SHEET
SKU: BTT145015 (Aluminium)
BTT145010 (Treated Steel)

Aspire Over Toilet Aid - Wide

PRODUCT OVERVIEW

- Raised armrests for easy transfers
- Height-adjustable legs with uneven surface adjustment
- Contoured, clip-on seat for comfort and easy cleaning
- Corrosion-resistant frame
- Removable splash guard for easy cleaning and transport
- Durable non-slip feet for stability in wet conditions

SPECIFICATIONS

Overall Height	700 - 880 mm	Overall Depth	625 mm
Overall Width	570 mm	Seat Depth	435 mm
Seat Width	370 mm	Seat Height (Adj.)	450 - 630 mm
Width Between Arms	525 mm	Unit Weight	3.4 kg
Max. Weight Capacity	175 kg	Material	Aluminium / Treated Steel