

Airgo Fusion Side-Fold Rollator and Transporter

SKU: 700-935

\$845



Features and Benefits

This Airgo Fusion is a two-in-one side-folding rollator & transport chair, to allow users independence and the ability to live life to its fullest. The Fusion's lightweight design and ergonomic-form handgrips provide stability and confidence.

Features

- Side-folding by the pull of a strap, for easy lifting and storage
- Padded seat and backrest for ultimate comfort
- Dual armrests for either mode
- Ergonomic handgrips provide easier brake access
- 20cm Comfort-Grip wheels are great for indoor and outdoor use
- Maximum user weight of 136kg

Specifications:

Product Width – 69cm

Product Depth – 72cm

Product Height - 76-93cm

Weight Capacity – 136kg

Handle Height Min – 76cm

Handle Height Max – 93cm

Seat to Floor Height – 51cm

Product Weight - 7.76kg

Seat Width – 50cm

Seat Depth – 28cm