

Ez Stand-N-Go

SKU: 43002

\$227



The leverage needed to help you stand

If you or a loved one has trouble standing from a couch, chair, or recliner, the EZ Stand & Go from Stander can help you stand with confidence; without having to buy an expensive lift chair. It's perfect for those who are recovering from surgery, have back pain, general unsteadiness, or have a lower-set couch that is hard to stand from.



Specifications:

Maximum User Weight:	136kg
Handle Height from Floor:	55.9cm - 81.3cm adjustable
Width Between Handles:	48.3cm - 66.0cm adjustable
Base Depth:	58.4cm - 66.0cm adjustable
Base Width:	48.3cm - 66.0cm adjustable
Tubing Diameter:	2.54cm
Product Weight:	3.9kg

