



# BME

BRIBIE MEDICAL EQUIPMENT

*Local Family  
Operated Business*

PRODUCT SHEET

SKU: KA410Z

## KCare Over Toilet Frame

### PRODUCT OVERVIEW

- Designed to assist users with limited mobility during toileting
- Adjustable height and angle to suit individual users and uneven floors
- Supportive arms to assist with sitting and transferring
- Clip-on plastic seat
- Removable splash guard
- Rubber leg tips provide anti-slip stability

### SPECIFICATIONS

|                |                                 |                      |              |
|----------------|---------------------------------|----------------------|--------------|
| Overall Height | 680 - 805 mm                    | Overall Depth        | 600 mm       |
| Overall Width  | 510 mm                          | Seat Height (Adj.)   | 480 - 605 mm |
| Seat Width     | 450 mm                          | Unit Weight          | 5 kg         |
| Material       | Epoxy Powder-Coated Steel Frame | Max. Weight Capacity | 110 kg       |