



BME

BRIBIE MEDICAL EQUIPMENT

*Local Family
Operated Business*

PRODUCT SHEET

SKU: RM402

RM Lightweight Aluminium Shower Chair



PRODUCT OVERVIEW

- Moulded plastic seat with drainage holes
- “S” bend side frame provides maximum space and comfort when seated
- Large rubber suction ferrule to avoid slipping and provide maximum safety

SPECIFICATIONS

Overall Height	790 - 940 mm	Overall Depth	640 mm
Overall Width	600 mm	Seat Depth	410 mm
Seat Width	490 mm	Seat Height (Adj.)	440 - 590 mm
Seat to Armrest Height	210 mm	Unit Weight	3.3 kg
Max. Weight Capacity	160 kg	Material	Aluminium