



BME

BRIBIE MEDICAL EQUIPMENT

*Local Family
Operated Business*

PRODUCT SHEET

SKU: RM400

RM Lightweight Over Toilet Aid

PRODUCT OVERVIEW

- Aluminium heavy-duty frame.
- Included splash guard.
- Privacy lid included.
- Easily removable seat and lid for cleaning.
- “S” bend side frame provides maximum space.
- Large rubber suction tips to prevent slipping.



SPECIFICATIONS

Overall Height	650 - 800 mm	Overall Depth	610 mm
Overall Width	580 mm	Seat Depth	430 mm
Seat Width	500 mm	Seat Height (Adj.)	440 - 590 mm
Armrest Height	210 mm	Unit Weight	3.6kg
Max. Weight Capacity	160 kg		