



BME

BRIBIE MEDICAL EQUIPMENT

Local Family

Operated Business

PRODUCT SHEET

SKU: BTS118005

Aspire Shower Chair

PRODUCT OVERVIEW

- Contoured seat with drainage holes for comfort and support
- Additional leg height adjustment for uneven surfaces
- Raised armrests to assist transfers
- Corrosion-resistant frame construction
- Durable non-slip feet for stability in wet conditions
- Enhanced seat drainage for hygiene and comfort

SPECIFICATIONS

Overall Height	735 - 910 mm	Overall Depth	570 mm
Overall Width	512 mm	Seat Depth	400 mm
Seat Width	440 mm	Seat Height (Adj.)	375 - 550 mm
Width between Arms	470 mm	Unit Weight	3.1 kg
Max. Weight Capacity	200 kg	Material	Aluminium