



BME

BRIBIE MEDICAL EQUIPMENT

*Local Family
Operated Business*

PRODUCT SHEET

SKU: BTT164050

Aspire Toilet Seat Raiser

With Swing Back Arms



PRODUCT OVERVIEW

- Raises toilet height to assist sitting and standing
- Available in 3 height options
- Universal fit for most standard toilets
- Durable, easy-clean moulded plastic design
- Swing-back arms for side transfers and easy access
- Adjustable mounting bolts for secure installation

SPECIFICATIONS

Overall Width	565 mm	Overall Height	370 mm
Seat Height	50 / 85 / 120 mm	Unit Weight	4.2 kg
Distance Between Arms	480 mm	Max. Weight Capacity	136 kg
Material	Moulded Plastic		