



BME

BRIBIE MEDICAL EQUIPMENT

*Local Family
Operated Business*

PRODUCT SHEET

SKU: RM403

RM Aluminium Shower Stool



PRODUCT OVERVIEW

- Moulded plastic seat with drainage holes
- “S” bend side frame provides maximum space and comfort when seated
- Large rubber suction ferrule to avoid slipping and provide maximum safety

SPECIFICATIONS

Overall Height	670 - 820 mm	Overall Depth	440 mm
Overall Width	580 mm	Seat Depth	330 mm
Seat Width	500 mm	Seat Height (Adj.)	490 - 640 mm
Seat to Armrest Height	210 mm	Unit Weight	2.8 kg
Max. Weight Capacity	160 kg	Material	Aluminium