



BME

BRIBIE MEDICAL EQUIPMENT

Local Family

Operated Business

PRODUCT SHEET

SKU: RG554H2

Redgum Space Saver Shower Stool

PRODUCT OVERVIEW

- Modern, functional, compact design
- Suitable for confined spaces
- Contoured seat with drainage holes
- Reversible handles with ribbed grips
- Durable, corrosion-resistant construction
- Easy to clean
- Flexible PVC leg tips

SPECIFICATIONS

Overall Height	640 - 765 mm	Overall Depth	450 mm
Overall Width	515 mm	Seat Depth	355 mm
Seat Width	470 mm	Seat Height (Adj.)	485 - 610 mm
Seat to Armrest Height	200 mm	Unit Weight	2.6 kg
Max. Weight Capacity	130 kg	Material	Aluminium