



BME

BRIBIE MEDICAL EQUIPMENT

*Local Family
Operated Business*

PRODUCT SHEET

SKU: BTS119005

Aspire Aluminium Shower Stool

PRODUCT OVERVIEW

- Easy clip-on seat
- Push-button height adjustment
- Adjustable legs for uneven surfaces
- Corrosion-resistant frame
- Raised armrests for transfers
- Anterior tilt for easier standing
- Comfortable seat for stiff or painful hips
- Non-slip feet for stability
- Suitable for wet conditions

SPECIFICATIONS

Overall Height	670 - 855 mm	Overall Depth	400 mm
Overall Width	545 mm	Seat Depth	310 mm
Seat Width	375 mm	Seat Height (Adj.)	470 - 655 mm
Width Between Arms	495 mm	Unit Weight	3.1 kg
Max. Weight Capacity	175 kg	Material	Aluminium