



BME

BRIBIE MEDICAL EQUIPMENT

*Local Family
Operated Business*

PRODUCT SHEET

SKU: BTS110900

Aspire Bath Transfer Bench - Padded

PRODUCT OVERVIEW

- Padded seat & backrest provides a cushioned surface for comfort
- Highly resilient vinyl covering
- Fully height adjustable with a side arm for secure transfers
- Adjustable backrest to suit right or left facing tubs
- Suction cup foot pieces secure the bench to the tub floor

SPECIFICATIONS

Seat Depth	405 mm	Safe Working Load	181 kg
Seat Width	600 mm	Material	Aluminium
Seat to Floor Height	455 - 555 mm		