



BME

BRIBIE MEDICAL EQUIPMENT

*Local Family
Operated Business*

PRODUCT SHEET

SKU: KA552

LEONARDO - Quad Motor Lift Recliner

PRODUCT OVERVIEW



- Dual-Lift Function
- Lay-Flat Recliner
- Footrest Extension
- All close & all lift buttons
- Ultra Durable Fabrics
- Lumbar support to accomodate sitting posture
- Four Motors
- Easy to use Hand Controller
 - soft grip buttons
 - 6 pre-defined positions
 - choose “full recline” or “full stand” with the touch of a button

SPECIFICATIONS

Overall Height	1092mm	Seat Height	483mm
Overall Width	902mm	Product Weight	76kg
Overall Depth	1067mm	Max User Weight	158kg
Seat Depth	521mm	Colour	Santa Barbara Fog
Seat Width	571mm		